


Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Pins hold a variety of wider-mouth flasks, beakers and other labware in inverted position. Retainer Top 4401901 is recommended for use with plastic and lightweight glassware (Top sold separately). Type 304 stainless steel pins have plastic tips to protect glassware from scratching. 19.3" w x 8.8" d x 8.0" h (49 x 22 x 20 cm).



Catalog Number: 4401801


Specifications

- **Weight:** 2.5 lbs
- **Weight metric:** 1.1 kg
- **Dimensions:** 19.3"w x 8.8"d x 6"h
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower, Upper

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

