

Overview



Pipette Inserts permit direct injection cleaning of pipets. Self-sealing silicone check valves allow washing of partial or full loads. Type 304 stainless steel.

6 Pipette Insert securely holds up to six pipettes vertically to clean and dry pipettes. Constructed of type 304 stainless steel. For use in Labconco Lower Spindle Racks manufactured after October 2020.

Catalog Number: 4575500


Specifications

- **Weight:** 3.0 lbs
- **Weight metric:** 1.4 kg
- **Dimensions:** 1.2" w x 8.7" d x 5.5" h
- **Dimensions metric:** 3.0 x 22.1 x 14.0 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

