


## Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Culture Tube Insert designed to hold 10-12 mm OD tubes. Retainer top included. Constructed of type 304 stainless steel. For use in Labconco Upper or Lower, Standard or Spindle Racks manufactured after October 2020.



**Catalog Number:** 4573600


## Specifications

- **Weight:** 5.0 lbs
- **Weight metric:** 2.3 kg
- **Dimensions:** 8.5" w x 18.5" d x 4.5" h
- **Dimensions metric:** 21.6 x 47.0 x 11.4 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower, Upper

## For Product Assistance



[Amy Ruano](#)  
Customer Service Representative



[Dan Simpson](#)  
Senior Sales Representative

