


Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Culture Tube Insert designed to hold 20-24 mm OD tubes. Constructed of type 304 stainless steel. For use in Labconco Lower Standard or Spindle Racks manufactured after October 2020.



Catalog Number: 4573800


Specifications

- **Weight:** 3.0 lbs
- **Weight metric:** 1.4 kg
- **Dimensions:** 8.5" w x 18.5" d x 9.8" h
- **Dimensions metric:** 21.6 x 47.0 x 24.9 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

