


Overview

Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Retainer Top secures culture tubes loaded into Culture Tube Inserts 4573600, 4573700, and 4573800. One is included with each Culture Tube Insert 10-12 mm size only (4573600). Constructed of type 304 stainless steel. For use in Labconco Upper or Lower, Standard or Spindle Racks manufactured after October 2020.



Catalog Number: 4573900


Specifications

- **Weight:** 1.5 lbs
- **Weight metric:** 0.7 kg
- **Dimensions:** 8.5" w x 19.1" d x 0.25" h
- **Dimensions metric:** 21.6 x 48.5 x 0.6 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower, Upper

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

