

Overview

Test Tube Holders are add-on accessories to the Lower Test Tube Spindle Rack. The Holders keep tubes upright. Constructed of Type 304 stainless steel.

Test Tube Spindle Conversion Kit includes six spindle arms with eight spindles per arm and a supportive base. Requires Test Tube Holder Insert 4574700, 4574800 or 4574900. For use in Labconco Lower Spindle Racks manufactured after October 2020.



Catalog Number: 4575300


Specifications

- **Weight:** 8.0 lbs
- **Weight metric:** 3.6 kg
- **Dimensions:** 7.2" w x 18.5" d x 2.3" h
- **Dimensions metric:** 18.3 x 47.0 x 5.8 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle
- **Rack Location:** Lower

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

