


Overview



Stainless steel Basket Inserts fit in Standard Open Racks and Spindle Racks.

Tray Insert for pans, instrument trays, multi-purpose trays and other shallow containers. Capacity: 5 trays. Constructed of stainless steel. For use in Labconco Lower Standard or Spindle Rack manufactured after October 2020.

Catalog Number: 4573000


Specifications

- **Weight:** 3.0 lbs
- **Weight metric:** 1.4 kg
- **Dimensions:** 9.5" w x 18" d x 9.3" h
- **Dimensions metric:** 24.1 x 45.7 x 23.6 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, FlaskScrubber Vantage Series Glassware Washers, SteamScrubber Glassware Washers
- **Applicable Rack Types:** Spindle, Standard
- **Rack Location:** Lower

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

