



Overview

Automatic dispenser kits and storage for detergent and neutralizing acid rinse with Labconco glassware washers.

Liquid Detergent Dispenser Kit connects to any SteamScrubber or FlaskScrubber washer to automate liquid detergent dispensing. Includes peristaltic pump, float switch, tubing, and one-gallon plastic container. Detergent not included. For use with washers manufactured after October 2020.

Catalog Number: 4587520


Specifications

- **Weight:** 3.0 lbs
- **Weight metric:** 1.4 kg
- **Dimensions:** 6.0" w x 6.0" d x 11.8" h
- **Dimensions metric:** 15.4 x 15.4 x 30.0 cm
- **Applicable Products:** FlaskScrubber Glassware Washers, SteamScrubber Glassware Washers

For Product Assistance



[Amy Ruano](#)
Customer Service Representative



[Dan Simpson](#)
Senior Sales Representative

